



TZ MEDICAL

RADIAL RUNWAY


Radial Runway

The Radial Runway provides a radiolucent platform that keeps the patient in normal anatomical position for ideal radial access. Disposable straps allow the wrist to be positioned in multiple neutral or extended angles optimizing patient comfort, and the foam pad may be removed from the base at the conclusion of the case. The disposable brace can then be attached to the foam pad, providing support for the patient's wrist during recovery.

Can be used on either side for radial access along with the TZ Medical Flat Radial Board. (Shown below.)



*Radial Runway
Reusable Base: Item# RRU-0008
Disposable Top Pad: Item# RRU-0009*



*The Radial Runway
in use with a Flat
Radial Board (Item#
CZ-400) for a radial
access procedure.*



Product Features

- Can be used for radial access on either wrist
- Adjustable straps fit any sized patient
- Disposable straps, foam pad, & brace
- Base can be sterilized and reused
- Keeps patient comfortable during procedure
- Detachable pad with brace remains with patient for recovery
- Latex Free

TZ Medical Inc.
17750 SW Upper Boones Ferry Rd. #150
Portland, OR 97224
Office: 800.944.0187 Fax: 503.639.0239
www.tzmedical.com

TZ Medical
Sparked by your ideas

